



KeyWords

ACCESS SYSTEM

SURGICAL STAPLING

SMOKE FILTRATION AND EVACUATION

RETRACTORS/STAYS & OTHER ACCESSORIES

HYDROGEL EYE PATCH

BREAST BIOPSY PROBE


DISPOSABLE ULTRASONIC SCALPEL

INTELLIGENT INSUFFLATOR SYSTEM

4K LAPAROSCOPY SYSTEM

BLADED TROCAR

HOME > PRODUCTS > ACCESS SYSTEM > BLADED TROCAR



Bladed Trocar Laparoscopic Disposable Trocars

Bladed Laparoscopic disposable trocar sharp linear blade allows low insertion force and reduced fascial defect.


Download

Send Inquiry

Bladed trocar Laparoscopic Disposable Trocars




reduced fascial defect.



5mm bladed trocar Structures



Appropriate Sealing offers good air tightness.



Threaded Design control insertion speed and depth.



Sharp linear blade allows low insertion force and reduced fascial defect.



Smooth conical aiming facilitates easier insertion.



Safety Lock indicates when the blade is shielded and exposed.



Simple seal detachment for versatility during specimen removal.

SURGAID Trocar laparoscopic surgery Specification

No.	Ref No.	Description
1	VC-100-1-5	bladed, 5*100mm
2	VC-100-1-10	bladed, 10*100mm
3	VC-100-1-12	bladed, 12*100mm
4	VCR-100-2-A	bladed, 2X5mm with smoke filter
5	VCR-100-2-B	bladed, 2X10mm with smoke filter
6	VCR-100-2-C	bladed, 2X12mm with smoke filter
7	VCR-100-3-A	bladed, 1X5mm, 2x10mm with smoke filter
8	VCR-100-4-A	bladed, 2X5mm, 2x10mm with smoke filter
9	VCR-100-4-B	bladed, 2X5mm, 1x10mm, 1x12mm with smoke filter
10	VCR-100-4-C	bladed, 2X5mm, 2x12mm with smoke filter

Note: For more product information, please consult our marketing department.

Trocars thoracic surgery Feature

1. One-time use to eliminate cross infection
2. Integrated design, no need to change the head, improve the operation efficiency;
3. Silicone valve, strong sealing, prevent air leakage;
4. Clear vision of transparent casing
5. The design of puncture head with new technology is less invasive and requires no wound suture
6. Three-stage sealed detachable design is more convenient to take things
7. Casing thread design has stronger adhesion

Applicable surgery types



Urology	Nephrectomy, adrenalectomy, ureterectomy, etc.
Gastrointestinal Surgery	Subtotal gastrectomy, appendectomy, ulcer perforation repair, intestinal adhesion lysis, colorectal tumor resection.
Thoracic Surgery	Lobectomy, esophageal cancer resection, etc.
General Surgery	Hernia repair, etc.


Certifications: CE, ISO



Shipment & Payment



WE SUPPORT THESE PAYMENT METHODS				
T / T	WESTERN UNION	L / C	PayPal	\$



air products. We

commit to seeing patients' benefits from the supply of minimally-invasive surgical devices, oriented towards a healthier life.

For Surgeons

Based on our surgeons' experience, we provide reliable and easy handling surgical trocar & instruments in order to improve the efficiency of the surgery and a better clinical experience.

For Distributors

Concentrate on market dynamics, Surgaid Medical provides distributors with competitive products and professional services. Thanks to our differential products and complete product line, our partners have an overwhelming advantage worldwide.

FAQ

Q: Can you a factory or a trading company?

A: We are a Chinese Medical instruments manufacturer.

Q: Can you accept the OEM?

A: Yes! we can do the OEM service for you.

Q: Which areas are you products usually sold to?

A: Our products usually exported to North America, South America, The Middle East, South-east Asia, Europe and so on.

Q: What' s your deliver time?

A: If we have the stock, it will be delivered with 3 days. If it is very large, it will be delivered within 8-45 days.

Hot Tags: bladed trocar laparoscopic disposable trocars, China, manufacturers, suppliers, factory, price,

Related Products




Intelligent 4k Laparoscopy System Camera



Laparoscopic Stapler Manufacturers



Disposable Linear Cutter Staplers



China Disposable Laparoscopic Trocar Set Factory



Disposable Single Access Endoscopic Surgery 4 Holes



Disposable Surgical Circular Circumcision Stapler

* Name:

* E-mail:



* Content: Please send us your enquiry details, we will reply you at the first time, thank you!

* Verification Code:

Send



Surgaid Medical (Xiamen) Co.,Ltd



Product Category

- Access System
- Surgical Stapling
- Accessory
- Smoke Filtration and Evacuation
- Hydrogel Eye Patch

CONTACT

- +86-592-6880461
- +86-592-6880677
- chenhui@surgaid.com
- No.218 Houxiang Road, Xinyang Industrial Zone, Haicang District, Xiamen, China

- Home
- About Us
- Products
- News
- Contact Us
- Shopping
- Mobile
- Sitemap

Copyright © Surgaid Medical (Xiamen) Co., Ltd. All Rights Reserved.Privacy Policy