



KeyWords

ACCESS SYSTEM

SURGICAL STAPLING

SMOKE FILTRATION AND EVACUATION

RETRACTORS/STAYS & OTHER ACCESSORIES

HYDROGEL EYE PATCH

BREAST BIOPSY PROBE

DISPOSABLE ULTRASONIC SCALPEL

INTELLIGENT INSUFFLATOR SYSTEM

4K LAPAROSCOPY SYSTEM

**BLADELESS TROCAR**

HOME > PRODUCTS > ACCESS SYSTEM > BLADELESS TROCAR



**Disposable Laparoscopic Insertion Single Use Trocar Kits**

bladeless disposable laparoscopic trocar in China is used to deliver CO<sub>2</sub> into abdominal cavity after abdominal paracentesis for thoracic laparoscopic surgery, and used in making channel for entry of endoscope and instruments from outside to abdominal cavity.

laparoscopic trocar and cannula

Send Inquiry

**Disposable Laparoscopic Insertion Single Use Trocar Kits**

**Product Description:**



nal paracentesis during thoracic  
n outside to abdominal cavity.We  
:ical single/four sheet(s) trocar  
rsal sealing design, large oriented  
f laparoscopic trocar china

manufacturer, please read down.  
disposable optical threaded design trocar :

### Disposable Laparoscopic Insertion Single Use Trocar Kits



**Universal Sealing design** increases resistance of instrument insertion without using the reducer.



**Large oriented aiming conical surface** facilitates one-handed operation and instrument insertion into the working channel.



**Threaded Design.** accepts both conventional puncture and spiral puncture. The controlled insertion speed and depth of trocar ensures stable fixation and avoids damage to abdomen.



**Translucent Design** for easier observation of instrument insertion/removal as well as extraction of specimen.




**Pyramid Tip** protects inter organs when puncturing.



**Additional Positioner** to ensure trocar placement .





#### Application:

Disposable Laparoscopic Trocar is extensively applied in various laparoscopic surgery including general minimally-Invasive surgery gynecological minimally-Invasive surgery, thoracic or urological endoscopic surgery, and more.

#### Disposable Laparoscopic Insertion Trocar Feature:

##### ✔Outstanding Ergonomics:

Ergonomic design offers increased comfort and better control  
Minimizes the risk of accidental injury

##### ✔Unique Protector:

It is effective to prevent from titanium clamp damaging the seal  
It can reduce the resistance of movement to the utmost extent

##### ✔Universal Seal:

Freely shift instruments between 12mm and 5mm with reliable closure



#### Laparoscopic Trocar Kits Specification:

No.	Ref No.	Description
1	NPCS-100-1-5	Bladeless, transparent cannula, 5*100mm
2	NPCS-100-1-10	Bladeless, transparent cannula, 10*100mm
3	NPCS-100-1-12	Bladeless, transparent cannula, 12*100mm
4	NPCS-100-1-15	Bladeless, transparent cannula, 15*100mm
5	NPCS-150-1-5	Bladeless, steel cannula, 5*150mm, For Bariatric Surgery
6	NPCS-150-1-10	Bladeless, steel cannula, 10*150mm, For Bariatric Surgery
7	NPCS-150-1-12	Bladeless, steel cannula, 12*150mm, For Bariatric Surgery

#### Why choose Surgaid Medical :

✈Continuous research and development of new products: For the development of international market, our company is to make two breakthroughs, namely network and service.



relationship with a number of hospitals around the world. We have made it to our top 10 partners. We are committed to quality management has won the favor



Stable quality: Depending on the needs of the user, Our company has been through professional quality system certification, CE certification and International Standard ISO 13485:2016, SFDA. Please feel free to purchase.

Q&A

Q: Can you do OEM and ODM?

A: Yes, OEM and ODM are both acceptable. The material, color, style can customize, the basic quantity we will advise after we discuss.

Q: Can we use our own logo?

A: Yes, we can print your private logo according to your request.

Q: Which certificate do you have?

A: CE and ISO13485.

Q: Which shipping way can you provide?

A: We can provide shipping by sea, by air and by express.

Q: What type of your product?


A: Surgaid offer a full range of high quality Surgical instruments such as Disposable Trocar, Single Port Laparoscopic, Smoke Filtration and Evacuation, Surgical Stapling and so on. including but not limited to the items listed above.

Hot Tags: disposable laparoscopic insertion single use trocar kits, China, manufacturers, suppliers, factory, price, Disposable Laparoscopic Trocar Set, Single Port Laparoscopic, Disposable Laparoscopic Trocar, Endoscopic Screw Trocar Set, Disposable Bariatric Trocar Instruments, 150 mm Bariatric Trocar

Related Products



Single Use Adult Male Foreskin Circumcision Clamps Kits




Disposable 12 mm Optical Trocar




Bladeless Disposable Bariatric Trocar Instruments



3mm Disposable Pediatric Trocar



Disposable Surgical Smoke Evacuation Equipment



Disposable Endo GIA Stapler-Green Reload

\* Name:

\* E-mail: Please input your email here, it is very important!

\* Phone/WhatsApp:



time, thank

you!

\* Verification Code:

Send



Surgaid Medical (Xiamen) Co.,Ltd



### Product Category

- Access System
- Surgical Stapling
- Accessory
- Smoke Filtration and Evacuation
- Hydrogel Eye Patch

### CONTACT

- +86-592-6880461
- +86-592-6880677
- chenhui@surgaid.com
- No.218 Houxiang Road, Xinyang Industrial Zone, Haicang District, Xiamen, China

Home About Us Products News Contact Us Share Us Mobile Sitemap

Copyright © Surgaid Medical (Xiamen) Co., Ltd. All Rights Reserved. Privacy Policy