



ACCESS SYSTEM

SURGICAL STAPLING

SMOKE FILTRATION AND EVACUATION

RETRACTORS/STAYS & OTHER ACCESSORIES

HYDROGEL EYE PATCH

BREAST BIOPSY PROBE

DISPOSABLE ULTRASONIC SCALPEL

INTELLIGENT INSUFFLATOR SYSTEM

4K LAPAROSCOPY SYSTEM

OPTICAL TROCAR

[HOME](#) > [PRODUCTS](#) > [ACCESS SYSTEM](#) > [OPTICAL TROCAR](#)



Disposable 5mm Optical Trocar

Disposable 5 mm Trocar Sleeves are available in high-quality stainless steel and high-performance plastic. With a working length of 100 mm, disposable 5mm optical trocar provides optimal access to the surgical site. Disposable 5mm Optical Trocar is CE marked and manufactured under ISO 13485 quality system. It is single-patient-use for puncturing abdominal cavity and supplying gas during laparoscopic surgery.

[Send Inquiry](#)

Laparoscopic 5mm trocar

Product Description:

Disposable 5 mm Trocar Sleeves are available in high-quality stainless steel and high-performance plastic. With a working length of 100 mm, disposable 5mm optical trocar provides optimal access to the surgical site. Disposable 5mm Optical Trocar is



turing abdominal cavity and



Appropriate Sealing offers good air tightness.



Threaded Design control insertion speed and depth.



Visualization tip allows excellent identification.



Smooth conical aiming facilitates easier insertion.



Simple seal detachment for versatility during specimen removal.



Fixed button with advanced fixation keep scope from moving during procedures.



Product Feature:

- Smooth conical aiming facilitates easier insertion
- Threaded design maximizes abdominal wall retention



No.	Ref No.	Description
1	NPVS-100-1-5	Optical, 5*100mm, steel cannula
2	NPVS-100-1-10	Optical, 10*100mm, steel cannula
3	NPVS-100-1-12	Optical, 12*100mm, transparent cannula
4	NPVS-100-2-B	Optical, Length-100mm (1 Trocar 10mm + 1 Cannula 10mm)
5	NPVS-100-3-A	Optical, Length-100mm (1 Trocar 5mm + 1 Trocar 10mm + 1 cannula 10mm)
6	NPVS-100-4-B	Optical, Length-100mm (1 Trocar 5mm + 1 Trocar 10mm + 1 Trocar 12mm + 1 Cannula 5mm)
7	NPVS-100-4-C	Optical, Length-100mm (1 Trocar 5mm + 1 Cannula 5mm + 1 Trocar 12mm + 1 cannula 12mm)

Basic info:

PRODUCT INFORMATION	
Product name	Disposable 5mm Optical Trocar
Brand	Surgaid medical
Application	Gynecology, Hepatobiliary Surgery, Urology, Gastrointestinal Surgery, Thoracic Surgery, General Surgery
Material:	Plastic
Feature	Single-use and optical
Origin	Xiamen, China
Type	Surgical trocar
Certificate	ISO13485,CE
Packages	blister+box+carton
MOQ	1 carton
Sterilize	EO Sterilize
Customized	OEM&ODM are welcome

Hot Tags: Disposable 5mm Optical Trocar, China, manufacturers, suppliers, factory, price,

Related Products



Single Use Adult Male Foreskin Circumcision Clamps Kits



Disposable Laparoscopic Insertion Single Use Trocar ...



Disposable Single Port Incision Laparoscopic Incisio...



Disposable Single Access Endoscopic Surgery 4 Holes

Disposable Circular PPH Surgical Staplers

Disposable Surgical Circular Circumcision Stapler

* Name:

* E-mail:

* Phone/WhatsApp:

Company:

Subject:

* Content:

* Verification Code:

Send



Surgaid Medical (Xiamen) Co.,Ltd



Product Category

- Access System
- Surgical Stapling
- Accessory
- Smoke Filtration and Evacuation
- Hydrogel Eye Patch

CONTACT

- +86-592-6880461
- +86-592-6880677
- chenhui@surgaid.com
- No.218 Houxiang Road, Xinyang Industrial Zone, Haicang District, Xiamen, China